

PROTOCOLO DE LAVADO DEL 1000D Y 600D

A continuación se detalla el protocolo de lavado para aquellas mochilas en 1000D y 600D que han sido manchadas con sangre u otros fluidos y haya que descontaminarlas:

- Se procede al vertido directo de peróxido de hidrógeno sobre las zonas manchadas.
- Se limpia con un papel antes de verter etanol al 70%.
- Aquellos productos sin estructura semirígida interna pueden ser lavados en lavadora doméstica a 40°C con un detergente para textil sin blanqueador óptico.
- Para una desinfección más profunda, se puede usar cualquier desinfectante apto para textiles.

1000D AND 600D WASHING PROTOCOL

Below is the washing protocol for those backpacks in 1000D and 600D that have been stained with blood or other fluids and have to be decontaminated:

- Proceed to apply hydrogen peroxide directly on the stained areas.
- Clean the area with a paper and apply 70% ethanol.
- Those products without semi-rigid inner structure, can be washed in a domestic washing machine at 40°C with a regular detergent without optical brightener.
- If further disinfection is needed, you can apply disinfectant product suitable for textile or textile-friendly.

ELITE BAGS, S.L.