

PROTOCOLO DE LAVADO DEL TARPAULIN

A continuación se detalla el protocolo de lavado para aquellas mochilas en tarpaulin que han sido manchadas con sangre u otros fluidos y haya que descontaminarlas:

- Lavar con agua y jabón neutro al 4%.
- Humedecer un paño con la solución de jabón y aplicar a la mochila para eliminar la suciedad.
- Secar con un papel.
- Para una desinfección más profunda, se puede usar cualquier desinfectante apto para textiles.

TARPAULIN WASHING PROTOCOL

Below is the washing protocol for those backpacks in tarpaulin that have been stained with blood or other fluids and have to be decontaminated:

- Wash using water and 4% neutral soap.
- Apply this solution on a cloth and wipe clean the backpack to remove dirt.
- Dry it with a paper towel or a clean cloth.
- If further disinfection is needed, you can apply textile or textile-friendly disinfectant.

ELITE BAGS, S.L.

Vial de los Cipreses, 9, 03006 Alicante, Spain.

Telf: (+34)902 212 210 · info@elitebags.es · www.elitebags.es